

Material Compatibility

14.10.2015

Product: **Descosept sensitive**
Mat. No.: **00-323DS-000**

The product has been tested according to SOP 04.02.012 and SOP 04.02.027 to be suitable with the following materials.

Recommended for the use on

Metals:

- Stainless Steel V2A
- Stainless Steel V4A
- Aluminum
- Brass
- Copper

Surfaces:

- Synthetic leather
- Polyurethan (PU)

Plastics:


- Polycarbonate (PC)
- Polyoxymethylene (POM)
- Acrylonitrile-butadiene-styrene (ABS)

■ As Recommended usage

■ As surfaces may vary in quality the product suitability should be checked by testing first on a small inconspicuous area.
To ensure materials compatibility please do not exceed recommended concentration, contact time and temperature.

Not recommended for the use on

■ Polymethyl methacrylate (PMMA) (Plexiglas/Altuglass)



Dr. Erik Brückner
Head of R&D Disinfection, **Dr. Schumacher GmbH**

